

Diffusion MRI Data Preprocessing
Eddy-current and Head Motion Correction
Align all raw images to the b0 image to correct for head motion and eddy current distortions
↓
EPI Induced Susceptibility Artifacts Correction
Correct the spatially nonlinear distortions caused by B ₀ inhomogeneities in Echo-planar imaging
↓
Brain Extraction
Remove voxels not necessary for analysis such as bone, dura, air, etc., leaving just the brain
↓
Reconstruct Local Diffusion Pattern
Fit a diffusion tensor model at each voxel on preprocessed and eddy current corrected data
↓
Tractography
Reconstruct brain connectivity graphs using brain tractography algorithms like FACT
↓
Brain Region Parcellation
Parcellate ROIs from T1-weighted structural MRI

Structural Brain Network Construction
Construct Network
Compute the network based on the generated label and the reconstructed whole brain tractography

FSL	AFNI	Free Surfer	Track Vis	3D Slider	Tortoise	MRtrix3	DSI Studio
✓	✓	✓		✓	✓	✓	✓
✓	✓	✓		✓	✓	✓	
✓	✓	✓		✓		✓	✓
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓		✓	✓		✓	✓
✓	✓				✓	✓	✓

Recommended Software: FSL, Metric, DSI Studio
--