

Diffusion MRI Data Preprocessing

Eddy-current and Head Motion Correction

Align all raw images to the b0 image to correct for head motion and eddy current distortions



EPI Induced Susceptibility Artifacts Correction

Correct the spatially nonlinear distortions caused by B_0 inhomogeneities in Echo-planar imaging



Brain Extraction

Remove voxels not necessary for analysis such as bone, dura, air, etc., leaving just the brain



Reconstruct Local Diffusion Pattern

Fit a diffusion tensor model at each voxel on preprocessed and eddy current corrected data



Tractography

Reconstruct brain connectivity graphs using whole brain tractography algorithms like FACT



Brain Region Parcellation

Parcellate ROIs from T1-weighted structural MRI and map those ROIs to DTI space

Structural Brain Network Construction

Construct Network

Compute the network based on the generated label and the reconstructed whole brain tractography

FSL	AFNI	Free Surfer	Track Vis	3D Slider	Tortoise	MRtrix3	DSI Studio	DIPY	Tracto Flow
✓	✓	✓		✓	✓	✓	✓	✓	✓
✓	✓	✓		✓	✓	✓			
✓	✓	✓		✓		✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
✓	✓		✓	✓		✓	✓	✓	✓
✓	✓				✓	✓	✓		✓

Recommended Software: FSL, Metric, DSI Studio