

Functional Neuroimaging Preprocessing

Brain Extraction (BE)

Remove voxels not necessary for analysis such as bone, dura, air, etc., leaving just the brain

Slice-Timing Correction (STC)

Adjust for the fact that each slice in the volume is taken at a different time, not all at once

Motion Correction/Realignment (MC/R)

Correct movement made during scanning by aligning all the functional images with one reference

Co-registration (CR)

Align participant's functional images with the anatomical structural images for localization

Normalization (NO)

Wrap the data across subjects to a template/atlas standardized space

Smoothing (SM)

Perform weighted averages of individual voxels with neighboring voxels

Functional Brain Network Construction

Brain Region Parcellation

Choose an anatomical atlas and segment each subject into the regions of interest defined by the given atlas

Network Construction

Given the ROI parcellation, calculate the edge by the time series activity correlation between different brain regions