

Functional MRI Data Preprocessing		SPM 12	AFNI	FSL	Free Surfer	CONN	fMRI Prep	ANTs	Nilearn
Brain Extraction		✓	✓	✓	✓		✓	✓	✓
Remove unnecessary voxels such as bone, air, etc. from T1, apply generated brain mask to fMRI data									
Slice-Timing Correction		✓	✓	✓	✓	✓	✓		
Adjust for the fact that each slice in the volume is taken at a different time, not all at once									
Motion Correction/Realignment		✓	✓	✓	✓	✓	✓	✓	
Correct movement made during scanning by aligning all the functional images with one reference									
Co-registration		✓	✓	✓	✓	✓	✓	✓	
Align participant’s functional images with the anatomical structural images for localization									
Normalization		✓	✓	✓	✓	✓	✓		
Warp the data across subjects to a template/atlas standardized space									
Smoothing		✓	✓	✓	✓	✓		✓	✓
Perform weighted averages of individual voxels with neighboring voxels									
Functional Brain Network Construction		Recommended Software: <b>CONN</b> , <b>GraphVar</b> , <b>Brain Connectivity Toolbox</b>							
Brain Region Parcellation	Construct Network								
Segment each subject into the ROI defined by the given atlas	Calculate pairwise correlations between ROIs as edges								