

## Question

Q263: What is the activity on the 100th day for an Ir-192 source ( $T_{1/2} = 78.4$  days) with an initial activity of 10 Ci?

- A. 2.8 Ci
- B. 4.1 Ci
- C. 5.1 Ci
- D. 7.8 Ci

## ChatGPT-3.5

To calculate the activity on the 100th day for an Ir-192 source, we can use the following formula:

$$A(100) = A(0) * 2^{(-t / T_{1/2})}$$

where  $A(0)$  is the initial activity,  $t$  is the time elapsed since the source was last calibrated, and  $T_{1/2}$  is the half-life of Ir-192.

Given that the half-life of Ir-192 is 78.4 days, the time elapsed since the source was last calibrated (i.e., the time between the calibration and the 100th day) is:

$$t = 100 \text{ days} - \text{calibration day}$$

We don't have the calibration day, so we cannot determine the exact value of  $t$ . However, we can calculate the fraction of the initial activity remaining on the 100th day as:

$$2^{(-t / T_{1/2})} = 2^{(-100 / 78.4)} = \mathbf{0.167}$$

Therefore, the activity on the 100th day can be calculated as:

$$A(100) = A(0) * 0.167$$

$$A(100) = 10 \text{ Ci} * 0.167$$

$$A(100) = 1.67 \text{ Ci}$$

Therefore, the activity on the 100th day for an Ir-192 source with an initial activity of 10 Ci is approximately 1.67 Ci.

Answer: None of the options are correct.

## ChatGPT-4

To find the activity of the Ir-192 source after 100 days, we can use the formula for radioactive decay:

$$\text{Activity\_final} = \text{Activity\_initial} * (1/2)^{(\text{time\_elapsed} / \text{half\_life})}$$

In this case, the initial activity ( $\text{Activity\_initial}$ ) is 10 Ci, the half-life ( $T_{1/2}$ ) is 78.4 days, and the time elapsed is 100 days. Plugging these values into the formula, we get:

$$\text{Activity\_final} = 10 * (1/2)^{(100 / 78.4)}$$

$$\text{Activity\_final} \approx 10 * (1/2)^{(1.276)}$$

$$\text{Activity\_final} \approx 10 * \mathbf{0.2835}$$

$$\text{Activity\_final} \approx 2.835 \text{ Ci}$$

The closest answer to this value is 2.8 Ci, so the correct choice is:

A. 2.8 Ci