Description of supporting movies

Movie S1: Periodic detachment of liquid lenses

Liquid: Dodecane Flow rate: $Q = 100 \,\mu\text{L/min}$

Injection of liquid 2 is done through a needle (its tip is visible at the top) resting against a vertical plate. The upper horizontal line is the triple contact line of the meniscus of liquid 1 (water) on the plate. The lower horizontal line is the image of the submerged bottom of the plate.

Movie S2: Example of periodic pinch-off (I) in the quasi-static limit i)

Liquid: Dodecane Flow rate: $Q = 5 \,\mu$ L/min Corresponding to Ca Bo₂ = 2.2 × 10⁻⁵

For quasi-static dripping (i), pinch-off occurs close below the contact line of the water meniscus on the vertical plate.

Movie S3: Example of periodic pinch-off (I) in the viscous limit ii)

Liquid: Mineral Oil S200 Flow rate: $Q = 50 \,\mu\text{L/min}$ Corresponding to Ca Bo₂ = 1.4×10^{-2}

For viscous pinch-off (ii), the floating liquid ligament stretches along the water meniscus into a thin thread, which eventually pinches-off.

Movie S4: Transition from periodic pinch-off (I) to no pinch-off (II)

Liquid: Mineral Oil RTM30 Flow rate: $Q = 20 \,\mu\text{L/min}$ Corresponding to Ca Bo₂ = 7.7 × 10⁻²

Liquid lenses detach, but the ligament width fluctuates before pinching-off. Pinch-off occurs when most of the lens volume rests on the horizontal part of the bath. The horizontal line is the image of the submerged bottom of the plate; the triple contact line (out of the video frame) is located slightly above.

Movie S5: No pinch-off regime (II) (continuous feeding)

Liquid: Mineral Oil RTM30 Flow rate: $Q = 50 \,\mu$ L/min Corresponding to Ca Bo₂ = 0.19

No pinch-off occurs, and a liquid ligament continuously feeds a floating liquid puddle. The horizontal line is the image of the submerged bottom of the plate; the triple contact line (out of the video frame) is located slightly above.